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* Messer, L., & Beatson, C., Preserving Skin Integrity with Chronic Device Use in Diabetes. Technology & Therapeutics Volume 20, Supplement 2, 2018.

† Kamann, S., Heinemann, L., & Oppel, E., Usage of Hydrocolloid-Based Plasters in Patients Who Have Developed Allergic Contact Dermatitis to Isobornyl Acrylate While Using Continuous Glucose Monitoring Systems. Journal of Diabetes Science and Technology, 2019.

‡ Paret, M., Barash, G. & Rachmiel, M. "Out of the box" solution for skin problems due to glucose-monitoring technology in youth with type 1 diabetes: real-life experience with fluticasone spray. Acta Diabetol 57, 419–424 (2020).

1. Ives B, Sikes K, Urban A, et al.: Practical aspects of realtime continuous glucose monitors: the experience of the Yale Children's Diabetes Program. Diabetes Educ 2010;36: 53–62.

2. Chase HP, Messer L: Understanding Insulin Pumps and Continuous Glucose Monitors. 3rd ed. Denver: Children's Diabetes Research Foundation, 2016.

3. Karlin AW, Ly TT, Pyle L, et al.: Duration of infusion set survival in lipohypertrophy versus nonlipohypertrophied tissue in patients with type 1 diabetes. Diabetes Technol Ther 2016;18:429–435.

4. McNichol L, Lund C, Rosen T, Gray M: Medical adhesives and patient safety: State of the science: consensus statements for the assessment, prevention, and treatment of adhesiverelated skin injuries. J Wound Ostomy Continence Nurs 2013;40:365–380; quiz E361–E362.

5. Englert K, Ruedy K, Coffey J, et al.: Skin and adhesive issues with continuous glucose monitors: a sticky situation. J Diabetes Sci Technol 2014;8:745–751.

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